



Mach7 EIP 2026.1

Enterprise Imaging Platform

IHE Integration Statement

Document ID: M7EIP-016

Revision date: July 30, 2026

Revision History

Document Version	Release Version	Release Date	Sections Affected	Description
2.4	2026.1	July 2026		Republished for Mach7 EIP v2026.1
2.3	12.5.0	April 2026		Republished for Mach7 EIP v12.5.0
2.2	12.4.0	August 2024		Republished for Mach7 EIP v12.4.0
2.1	12.3.0	June 2024		Republished for Mach7 EIP v12.3.0
2.0	12.2.0	February 2024	All	Document refreshed, new template and formatting applied. Removed Cross-Enterprise Document Sharing (XDS) integration profiles.
1.1	12.2.0	November 2023	General Release	General release for Mach7 EIP v12.2.0
1.0	12.0.0	November 2022	General Release	General release for Mach7 EIP v12.0.0

IHE Technical Framework Integration

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles and Actors listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Audit Trail and Node Authentication (ATNA)	Secure Node	none
Access to Radiology Information (ARI)	Image Manager	none
	Image Archive	none
	Report Repository	none
Consistent Presentation of Images (CPI)	Image Manager	none
	Image Archive	none
Consistent Time (CT)	Time Client	none
Evidence Documents (ED)	Image Manager	none
	Image Archive	none
Image Object Change Management (IOCM)	Image Manager	none
	Image Archive	none
	Image Display	none
	Change Requester	none
Key Image Note (KIN)	Image Manager	none
	Image Archive	none

Integration Profiles Implemented	Actors Implemented	Options Implemented
Mammography Image (MAMMO)	Image Manager	none
	Image Archive	none
Multiple Image Manager Archive (MIMA)*	Image Manager	none
	Image Archive	none
Nuclear Medicine Image (NMI)	Image Manager	none
	Image Archive	none
Patient Identifier Cross-Reference HL7 V3 (PIXV3)	Identity Consumer	none
Patient Information Reconciliation (PIR)*	Image Manager	none
	Image Archive	none
Portable Data for Imaging (PDI)	Image Display	none
	Portable Media Importer	none
	Portable Media Creator	none
Radiation Exposure Monitoring (REM)	Image Manager	none
	Image Archive	none
Scheduled Workflow (SWF)*	Image Manager	none
	Image Archive	none
Teaching File and Clinical Trials (TCE)	Export Manager	none

Link to Standards Conformance Statements for the Implementation

HL7	https://www.mach7t.com/collateral
DICOM	