

eUnity

by  MACH7

eUnity 7.10

Enterprise Viewing, Integration, and Collaboration Platform

IHE Integration Statement

Document ID: EUNITY-024

Revision number: 1.8

Revision date: July 24, 2026

Table of contents

Revision history	3
IHE integration statement	4
Links to general information on IHE	4
eUnity DICOM conformance	4

Revision history

Revision number	Date	Description of changes
1.8	July 24, 2026	Updated for 7.10
1.7	May 15, 2026	Updated for 7.9 Updated version number.
1.6	November 19, 2025	Updated for 7.8 Updated version number.
1.5	May 7, 2025	Updated for 7.7. Updated version number.
1.4	February 5, 2025	Updated for 7.6. Updated version number.
1.3	June 17, 2024	Updated for 7.5. Updated version number.
1.2	October 3, 2023	Updated for 7.4. Updated version number.
1.1	March 22, 2023	Updated for 7.3. Updated version number.
1.0	July 7, 2022	Updated for 7.2. Updated version number.

IHE integration statement

Product: eUnity

Vendor: Mach7 Technologies Canada Inc.

Integration profiles implemented	Actors implemented	Options implemented
Scheduled Workflow (SWF)	Image display	None
Consistent Presentation of Images (CPI)	Image display	None
Simple Image and Numeric Report (SINR)	Report reader	None
Access to Radiology Information (ARI)	Report reader	Multiple sources
	Image display	Multiple sources
Audit Trail and Node Authentication (ATNA)	Secure application	Encryption option
	Audit record repository	None

Links to general information on IHE

North America	www.ihe.net
Europe	www.ihe-europe.net
Japan	http://www.ihe-j.org/

eUnity DICOM conformance

For the eUnity DICOM conformance statement, see <https://mach7t.com/collateral>.